









## Cryopreparation – Stem Cell Banking

### Cryopreparation package

Biosafe offers a wide range of cryopreparation accessories to ensure a safe and efficient preparation of the cryobag containing concentrated hematopoietic stem cells before freezing. With these accessories the complete workflow of cryopreparation is covered.

Products included:			
		Traceability labels for cryopreservation or Traceability labels for cryopreservation – BioArchive® System format	4-label set comprising of label(s) for the cryobag (single or double compartment), for the canister and closing seal (safety closure)
		Overwrap bag	Overwrap to protect the cryobag during freezing process, made by EVA material, 95 x 92 mm (3.7" x 3.6"). It remains flexible in liquid nitrogen.  To be used with Biosafe overwrap sealing device
		Insertion tool	Used to facilitate the insertion of the cryobag into the overwrap bag
		Metal storage canister	Aluminium alloy canister ensures enhanced thermal conductivity for homogeneous freezing. External dimensions: 94.4x80.4x10.5 mm (3.7"x3.2"x0.4"). Internal: 91.5 x 77.0 x 9.6 mm (3.6"x3.0"x0.4")
		AK-100 Sterile Sampling Line	Sampling line for product cell-count. AK-100 is connected to the product bag using SCD device. For samples up to 3 ml
		AK-101 Sterile Sampling Line	Sampling line for product cell-count. AK-101 is connected to the product bag using SCD device. For samples up to 8 ml
		FB-100.1/2b Single/double compartment Cryobags	Stand-alone, sterile cryobags. Volume capacity of 10-30 ml for CryoSC-S and 25 ml for CryoSC-Db; connection by spike or SCD, and port for DMSO injection
		FA-100.1 DMSO extension line	To allow extended line for DMSO addition

## Order information

Traceability labels for cryopreservation

Product #7011 (Biosafe format) or #7031 (BioArchive<sup>®</sup> format)  
Minimum order quantity is 1 roll of 1000 sets of labels

Overwrap bag

Product #4085, minimum order quantity is 1 box of 100 units

Insertion tool

Product #7030, minimum order quantity is 1 package of 2 units

Metal storage canister

Product #7012, minimum order quantity is 1 box of 40 units

AK-100 Sterile Sampling Line

Product #10010, minimum order quantity is 1 box of 60 units

AK-101 Sterile Sampling Line

Product #10028, minimum order quantity is 1 box of 60 units

Stand-alone CryoSC bag

Product #10023 for CryoSC-S or #10041 for FB-100.2b

Minimum order quantity is 1 box of 24 units

DMSO extension line

Product #10022, minimum order quantity is 1 box of 24 units

All cryopreparation accessories are CE marked.

Head Office  
**Biosafe Group SA**

Eysins, Switzerland  
T: +41 22 365 27 27  
E: info@biosafe.ch

North America Office  
**Biosafe America, Inc.**

Houston, USA  
T: +1 713 936 2615  
E: info@biosafeamerica.com

Latin America office  
**Biosafe Latin America**

São Paulo, Brazil  
T: +55 11 2655 1731  
E: info@biosafelatinamerica.com

Asia office  
**Biosafe Asia-Pacific Ltd.**

Hong Kong  
T: +852 2956 7500  
E: info@biosafeasiapacific.com

China office  
**Biosafe Medical Device Int. Trading  
(Shanghai) Co. Ltd**

Shanghai, China  
T: +86 21 6137 3214  
E: info@biosafechina.com