



新揚行+

Cosmo Medical Supply

高雄市左營區明誠二路332號4樓之3

TEL: (07)5582188 FAX: (07)5569943

<https://www.cosmomed.com.tw>



VASINI STRUMENTI S.r.l.

Divisione Elettromedicale

# AUTO SEALER


REF : AS-03

PORTABLE SEALER **FOR**  
**PVC and EVA tubes**

**More than 2000 sealing cycles WITH FULL BATTERIE**



Battery charge  
indicator



OPTION  
Elastomer ring  
avoids sliding



Allow sealing close  
to the blood bag



Ergonomic design for  
easy cleaning



Visual LED  
displays sealing cycle

- Easy to use
- Constant sealing
- Automatic recharge of battery less than 2 hours
- Low power consumption
- Light weight (less than 1,25 kg)
- Constant sealing quality
- Level of battery charge controlled by 6 LED
- Able to be used during the blood donation process




High quality of sealing and easy separation

+390544 50 161 7

[export@vasinistrumenti.it](mailto:export@vasinistrumenti.it)


VASINI STRUMENTI

- Total control of sealing process by Micro-Processor
- Ecosystem friendly "Lithium Ion Phosphate" battery and high capacity for safety and extended lifetime



Seal all kind of PV blood bag tubes (min 3mm and up to 6 mm, with thickness of up to 0,8 mm) with motorized clamping system

## TECHNICAL FEATURES

<b>Main power:</b> AC/DC charger <b>Input:</b> 100-240 VAC, 47-63 Hz, 0,3 - 0,9A Max <b>Output:</b> 12 Vdc / 5,25A max <b>Battery:</b> Rechargeable battery 14,8 V 3000 mA Lithium-ION latest technology <b>RF output max:</b> 90 W - 50 ohm <b>RF frequency:</b> 40,685 MHz <b>Operating temperature:</b> 10°C to 45°C (50°F à 113°F) <b>Humidity:</b> 10% - 90% (without condensation) <b>Storage temperature:</b> -10°C to 50°		
<b>BASE (Aluminium)</b> <b>Dimensions :</b> 23 x 8,5 x 4 cm <b>Weight</b> 0,75 kg	<b>CARRYING CASE</b> <b>Dimensions :</b> 44,5 x 37 x 11 cm <b>Weight:</b> 1,20 kg	
<b>CLAMP (Aluminium)</b> <b>Dimensions :</b> 24 cm (long), Ø 3,5 mm <b>Weight</b> 0,33 kg	<b>PACKING</b> <b>Dimensions:</b> 47 x 39 x 14 cm <b>Weight:</b> 2,54 kg (with case) / 1,23 kg (without case)	
<b>Sealing time:</b> 2 seconds		
<b>CERTIFICATIONS:</b> ISO 9001, Calibration and electrical certificate		<b>CLASS I TYPE B</b>

## ACCESORIES

		
Transport suitcase	Suitcase for 3 units	Battery

### Blood Bank Supplies

Via Faentina n. 207/Q – 48124 RAVENNA – ITALY T: +39  
0544 46404 F: +39 0544 501617

[www.vasinistrumenti.it](http://www.vasinistrumenti.it) [export@vasinistrumenti.it](mailto:export@vasinistrumenti.it)

